

Supplementary Table S1 Comparison data between indirect immunofluorescence (reference method) and ELISA (comparative method) with binary results (positive/negative) in detection of ANCA

ELISA		Indirect immunofluorescence				
		Overall positive	Overall negative	C-ANCA	P-ANCA (FS)	P-ANCA (FR)
Overall positive		20	2	10	6	4
Overall negative		13	154	10	0	3
PR3-ELISA	Positive ^a	10	1	8	2	0
	Negative	23	155	12	4	7
MPO-ELISA	Positive ^a	11	1	2	3	6
	Negative	22	155	18	3	1

Abbreviations: ANCA, antineutrophil cytoplasmic antibody; C-ANCA, cytoplasmic pattern of ANCA; ELISA, enzyme-linked immunosorbent assay; FR, formalin resistant; FS, formalin sensitive; MPO, myeloperoxidase; P-ANCA, perinuclear pattern of ANCA; PR3, proteinase 3.

^aPositive indicates ≥ 20 units.